




SILICA WARE

1. LABORATORY SILICA WARE

1.1 Round Bottom Basins with Spout

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
S-600-100	Round Bottom Basins, 20	23	55
S-600-200	Round Bottom Basins, 45	27	75
S-600-300	Round Bottom Basins, 70	33	87
S-600-400	Round Bottom Basins, 85	36	94
S-600-500	Round Bottom Basins, 100	38	100
S-600-600	Round Bottom Basins, 200	47	113
S-600-700	Round Bottom Basins, 400	62	142



1.2 Deep Form Basin with Spout

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
S-600-800	Deep Form Basin, 10	18	40
S-600-900	Deep Form Basin, 15	25	50
S-601-000	Deep Form Basin, 45	30	60
S-601-100	Deep Form Basin, 90	45	80
S-601-200	Deep Form Basin, 170	55	100



1.3 Flat Form Basins with Spout

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
S-601-300	Flat Form Basins with Spout 20 ml	15	69
S-601-400	Flat Form Basins with Spout 25ml	15	79
S-601-500	Flat Form Basins with Spout 75ml	20	101
S-601-600	Flat Form Basins with Spout 125ml	23	131



1.4 Low Form Crucibles and Lids

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
S-601-700	Silica Crucibles 5ml Without lid	19	29
L-601-700	Lid For Cat. No. S-601-700	-	-
S-601-800	Silica Crucibles 10ml Without lid	19	38
L-601-800	Lid For Cat. No. S-601-800	-	-
S-601-900	Silica Crucibles 15ml Without lid	25	41
L-601-900	Lid For Cat. No. S-601-900	-	-
S-602-000	Silica Crucibles 25ml Without lid	28	47
L-602-000	Lid For Cat. No. S-602-000	-	-
S-602-100	Silica Crucibles 50ml Without lid	37	57
L-602-100	Lid For Cat. No. S-602-100	-	-
S-602-200	Silica Crucibles 80ml Without lid	45	67
L-602-200	Lid For Cat. No. S-602-200	-	-
S-602-300	Silica Crucibles 150ml Without lid	51	82
L-602-300	Lid For Cat. No. S-602-300	-	-
S-602-400	Silica Crucibles 100ml Without lid	48	72
L-602-400	Lid For Cat. No. S-602-400	-	-
S-602-500	Silica Crucibles 125ml Without lid	-	-
L-602-500	Lid For Cat. No. S-602-500	-	-
S-602-600	Silica Crucibles 250ml Without lid	-	-
L-602-600	Lid For Cat. No. S-602-600	-	-



1.5 Tall Form Crucibles & Lids

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
S-602-700	Silica Crucibles 10ml Without lid	32	27
L-602-700	Lid For Cat. No. S-602-700	-	-
S-602-800	Silica Crucibles 20ml Without lid	38	35
L-602-800	Lid For Cat. No. S-602-800	-	-
S-602-900	Silica Crucibles 30ml Without lid	39	43
L-602-900	Lid For Cat. No. S-602-900	-	-
S-603-000	Silica Crucibles 50ml Without lid	51	51
L-603-000	Lid For Cat. No. S-603-000	-	-
S-603-100	Silica Crucibles 200ml Without lid	81	73
L-603-100	Lid For Cat. No. S-603-100	-	-



1.6 Broad Base Crucibles

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	External Dimensions		
	Nominal Capacity (ml)	Height H (mm)	Diameter D (mm)
S-603-200	Broad Base Crucibles, 15	24	41
S-603-300	Broad Base Crucibles, 25	27	48
S-603-400	Broad Base Crucibles, 40	37	57



1.7 Boats Combustion with Handle

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	External Dimensions		
	Length L (mm)	Width W (mm)	Height H (mm)
S-603-500	Boats Combustion with Handle, 49	15	9
S-603-600	Boats Combustion with Handle, 77	12	10
S-603-700	Boats Combustion with Handle, 77	17	11
S-603-800	Boats Combustion with Handle, 77	20	12
S-603-900	Boats Combustion with Handle, 103	20	12



1.8 Moisture & Ash Determination Circular Capsule with Ground Lid (B.S 705)

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	External Dimensions	
	Height H (mm)	Diameter D(mm)
S-604-000	11	38
S-604-100	12	54
S-604-200	12	75



1.9 Volatile Matter Determination Coking Crucible & Lid B.S. 1016 (Part 1)

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Description	External Dimensions	
		Width W (mm)	Height H (mm)
S-604-300	Crucible with projection, Without Lid	38	25
S-604-400	Crucible without projection, Without Lid	38	25
S-604-500	Lids Use for S-604-300 & S-604-400	3 - 4.75	Di-21 to suit the crucible
S-604-600	Plunger Use for S-604-300 & S-604-400	12.5	to suit the crucible



1.10 Muffle Trays

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Length L (mm)	No. of holes
S-604-700	Muffle Trays 39 mm	2
S-604-800	Muffle Trays 76 mm	4



1.11 Capsule Circular

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
S-604-900	Capsule Circular 10 ml	15	40
S-605-000	Capsule Circular 20 ml	15	48
S-605-100	Capsule Circular 25 ml	15	54
S-605-200	Capsule Circular 30 ml	15	62
S-605-300	Capsule Circular 35 ml	15	69
S-605-400	Capsule Circular 50 ml	18	78
S-605-500	Capsule Circular 45 ml	28	56
S-605-600	Capsule Circular 50 ml	28	62
S-605-700	Capsule Circular 80 ml	28	75
S-605-800	Capsule Circular 100 ml	28	87



1.12 Tobacco Ash Capsules

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions		
		Height H (mm)	Width W (mm)	Length L (mm)
S-605-900	Tobacco Ash Capsules 90 ml	28	56	78



1.13 Rectangular Capsules

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions		
		Height H (mm)	Width W (mm)	Length L (mm)
S-606-000	Rectangular Capsules, 10	12	28	60
S-606-100	Rectangular Capsules, 20	12	39	68
S-606-200	Rectangular Capsules, 25	18	42	52



1.14 Gooch Crucible, Translucent

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Height H (mm)	Diameter D (mm)
S-606-300	Gooch Crucible, 52	60



1.15 Triangle on Nichrome Wire

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Length of side L (mm)	Nominal Diameter of inscribed circle D (mm)
S-606-400	Triangle on Nichrome Wire 38 mm	20 mm
S-606-500	Triangle on Nichrome Wire 44 mm	24 mm
S-606-600	Triangle on Nichrome Wire 50 mm	28 mm
S-606-700	Triangle on Nichrome Wire 58 mm	32 mm
S-606-800	Triangle on Nichrome Wire 63 mm	35 mm
S-606-900	Triangle on Nichrome Wire 69 mm	38 mm
S-607-000	Triangle on Nichrome Wire 76 mm	43 mm



2. INDUSTRIAL SILICA WARE

2.1 Basins, Round Bottom with Spout

- Made of silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
IS-607-100	Basins Round Bottom with Spout, 400	62	140
IS-607-200	Basins Round Bottom with Spout, 500	66	140
IS-607-300	Basins Round Bottom with Spout, 600	70	160
IS-607-400	Basins Round Bottom with Spout, 800	76	185
IS-607-500	Basins Round Bottom with Spout, 1200	100	185



2.2 Basins, Round Bottom without Spout

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
IS-607-600	Basins, Round Bottom without Spout, 1800	92	190



2.3 Basins, Flat Bottom with Spout

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
IS-607-700	Basins, Flat Bottom with Spout, 700	75	165



2.4 Crucibles, Tall Form without Spout

- Made of silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
IS-607-800	Crucibles, Tall Form without Spout, 900	155	115
IS-607-900	Crucibles, Tall Form without Spout, 2500	213	159



2.5 Crucibles Tall Form with Spout

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Overall Height H (mm)	External Diameter D (mm)
IS-608-000	Crucibles Tall Form with Spout, 350	111	95
IS-608-100	Crucibles Tall Form with Spout, 500	120	108
IS-608-200	Crucibles Tall Form with Spout, 900	155	115
IS-608-300	Crucibles Tall Form with Spout, 1500	190	125
IS-608-400	Crucibles Tall Form with Spout, 2500	213	159



2.6 Crucibles Conical Form without Spout

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Nominal Capacity (ml)	External Dimensions	
		Height H (mm)	Diameter D (mm)
IS-608-500	Crucibles Conical Form without Spout, 2700	305	150



2.7 Trays, Rectangular Trays

- Made of Silica




Ordering Information with Specifications:

Catalogue No.	Length L (mm)	Breadth B (mm)	Thickness T (mm)
IS-608-600	Trays Rectangular Trays, 75	75	25
IS-608-700	Trays Rectangular Trays, 100	75	25
IS-608-800	Trays Rectangular Trays, 150	75	25
IS-608-900	Trays Rectangular Trays, 100	100	25
IS-609-000	Trays Rectangular Trays, 150	100	25
IS-609-100	Trays Rectangular Trays, 100	100	30
IS-609-200	Trays Rectangular Trays, 150	100	30
IS-609-300	Trays Rectangular Trays, 150	150	25
IS-609-400	Trays Rectangular Trays, 150	150	30



2.8 Plates, Translucent Silica Plate

- Made of Silica



Ordering Information with Specifications:

Catalogue No.	Length L (mm)	Breadth M (mm)	Thickness T (mm)
IS-609-500	Plates, Translucent Silica Plate, 75	75	2mm
IS-609-600	Plates, Translucent Silica Plate, 75	75	3mm
IS-609-700	Plates, Translucent Silica Plate, 75	75	5mm
IS-609-800	Plates, Translucent Silica Plate, 100	100	2mm
IS-609-900	Plates, Translucent Silica Plate, 100	100	3mm
IS-610-000	Plates, Translucent Silica Plate, 100	100	5mm
IS-610-100	Plates, Translucent Silica Plate, 150	150	2mm
IS-610-200	Plates, Translucent Silica Plate, 150	150	3mm
IS-610-300	Plates, Translucent Silica Plate, 150	150	5mm