



SHANGHAI RONGTAI BIOCHEMICAL ENGINEERING CO., LTD.

Website: www.rt-bio.com Email: sales@rt-bio.com

Address: Building 8, No 255, Cao Nong Road, Song Jiang District, Shanghai, 201612, P. R. China

Tel: 8621 3701 3385 · 3701 3386

Fax: 8621 5988 2370



Product Catalog

SHANGHAI RONG TAI BIOCHEMICAL ENGINEERING CO., LTD.




About us

Shanghai Rong Tai Biochemical Engineering Co., Ltd. is a company which is specialized in the manufacture and R&D of lab analysis instruments and lab consumptive material.

We have product of Pipette, Bottle-Top Dispense, Centrifuge and microliter syringe etc. With more than 20 years professional experience we have got the license for making calibration device which issued by State Quality Supervision Bureau and have past International quality management control system, acquired the ISO9001 and CE. We have a complete range of specifications, provide comprehensive and extensive choice for lab research and clinical diagnosis.

We are always insist on the principle of quality first, customer for most. We gain the trust from our customers by stable quality and perfect service system. Our service objective is devoted to provide high added value product and professional service for all customers and partners.



Contents

Single Channel Digital Variable Pipette	03
Multi-channel Pipette	05
Pipetfocus	
Whole Autoclaved Single Channel Adjustable Pipette	07
Pipette with switch	09
Coloured Adjustable Pipette-Five Fixed Volume	11
Purple Adjustable Pipette-Five Fixed Volume	13
Dispensing Pipette	15
Adjustable Glass-Injection Dispenser	17
Bottle-Top Dispenser	18
Pen PH Recording	19
Microliter Syringe	20
Enzyme Label Plate	21
Tips for Pipette	22
Racks for Tips & Centrifuge Tubes	23
Pipette Stander& Tube Rack	24
Centrifuge Tubes	25
Others	26

Single Channel Digital Variable Pipette



An ergonomic design, convenience of use and simplicity of operation
 Ease of calibration and maintenance
 Availability for autoclave of part of the device
 Detailed instructions in Chinese and English
 A hook, disassembly tool and consumable replacement parts

- An ergonomic handle hook design for easy grip in the palm
- A clever design for the tip eject button
- A clear reading window
- Availability for autoclave of the tip and the pusher, also detachable to facilitate cleaning and maintenance
- A tip cone made of durable chemical materials, with excellent chemical corrosion resistance.



The Parameters of Single Channel Digital Variable Pipette

Product Code	Channel	Volume Range	Increment	Test Volume	Accuracy	Precision	Tips Code
20001010	1	0.5-10ul	0.1ul	10ul	±1.00%	±0.80%	80102010
				5ul	±1.50%	±1.50%	
				1ul	±2.50%	±1.50%	
20001015	1	2-20ul	0.5ul	20ul	±0.90%	±0.40%	80102060
				10ul	±1.20%	±1.00%	
				2ul	±3.00%	±2.00%	
20001020	1	5-50ul	0.5ul	50ul	±0.60%	±0.30%	80102060
				25ul	±0.90%	±0.60%	
				5ul	±2.00%	±2.00%	
20001025	1	10-100ul	1ul	100ul	±0.80%	±0.15%	80102060
				50ul	±1.00%	±0.40%	
				10ul	±3.00%	±1.50%	
20001030	1	20-200ul	1ul	200ul	±0.60%	±0.15%	80102060
				100ul	±0.80%	±0.30%	
				20ul	±3.00%	±1.00%	
20001040	1	100-1000ul	5ul	1000ul	±0.60%	±0.20%	80102090
				500ul	±0.70%	±0.25%	
				100ul	±1.50%	±0.70%	
20001050	1	1000-5000ul	50ul	5000ul	±0.50%	±0.15%	80102100
				2500ul	±0.60%	±0.30%	
				1000ul	±0.70%	±0.30%	

Order Details of Single Channel Digital Variable Pipette

Product Code	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
20001010	Single Channel Digital Variable Pipette	0.5-10ul	1	pcs	68.00	50pcs/carton	9.5	0.475*0.39*0.28
20001015	Single Channel Digital Variable Pipette	2-20ul	1	pcs	68.00	50pcs/carton	9.5	0.475*0.39*0.28
20001020	Single Channel Digital Variable Pipette	5-50ul	1	pcs	68.00	50pcs/carton	9.5	0.475*0.39*0.28
20001025	Single Channel Digital Variable Pipette	10-100ul	1	pcs	68.00	50pcs/carton	9.7	0.475*0.39*0.28
20001030	Single Channel Digital Variable Pipette	20-200ul	1	pcs	68.00	50pcs/carton	9.75	0.475*0.39*0.28
20001040	Single Channel Digital Variable Pipette	100-1000ul	1	pcs	68.00	50pcs/carton	10.2	0.475*0.39*0.28
20001050	Single Channel Digital Variable Pipette	1000-5000ul	1	pcs	68.00	50pcs/carton	10.3	0.475*0.39*0.28
20002000	Single Channel Digital Fixed Pipette	0.5ul	1	pcs	50.00	50pcs/carton	9.5	0.475*0.39*0.28
20002005	Single Channel Digital Fixed Pipette	1ul	1	pcs	50.00	50pcs/carton	9.5	0.475*0.39*0.28
20002010	Single Channel Digital Fixed Pipette	5ul	1	pcs	50.00	50pcs/carton	9.5	0.475*0.39*0.28
20002020	Single Channel Digital Fixed Pipette	10ul	1	pcs	50.00	50pcs/carton	9.5	0.475*0.39*0.28
20002030	Single Channel Digital Fixed Pipette	20ul	1	pcs	50.00	50pcs/carton	9.5	0.475*0.39*0.28
20002040	Single Channel Digital Fixed Pipette	25ul	1	pcs	50.00	50pcs/carton	9.75	0.475*0.39*0.28
20002050	Single Channel Digital Fixed Pipette	50ul	1	pcs	50.00	50pcs/carton	9.75	0.475*0.39*0.28
20002060	Single Channel Digital Fixed Pipette	100ul	1	pcs	50.00	50pcs/carton	9.75	0.475*0.39*0.28
20002070	Single Channel Digital Fixed Pipette	200ul	1	pcs	50.00	50pcs/carton	9.75	0.475*0.39*0.28
20002080	Single Channel Digital Fixed Pipette	250ul	1	pcs	50.00	50pcs/carton	10.2	0.475*0.39*0.28
20002090	Single Channel Digital Fixed Pipette	500ul	1	pcs	50.00	50pcs/carton	10.2	0.475*0.39*0.28
20002100	Single Channel Digital Fixed Pipette	1000ul	1	pcs	50.00	50pcs/carton	10.2	0.475*0.39*0.28
20002110	Single Channel Digital Fixed Pipette	2000ul	1	pcs	50.00	50pcs/carton	10.3	0.475*0.39*0.28
20002120	Single Channel Digital Fixed Pipette	2500ul	1	pcs	50.00	50pcs/carton	10.3	0.475*0.39*0.28
20002130	Single Channel Digital Fixed Pipette	5000u	1	pcs	50.00	50pcs/carton	10.3	0.475*0.39*0.28

Multi-channel Pipette



8 channel pipettes are available for standard 96-well plate
 Dispensing head rotates for optimum pipetting convenience
 Individual piston and tip cone assemblies allow easy repair and maintenance
 Compound material tip cone design allows visual seal verification
 Can be used with universal style pipette tips



The Parameters of Multi-channel Pipette

Product Code	Channel	Volume Range	Increment	Test Volume	Accuracy	Precision	Tips Code
20001110	8-channel	0.5-10ul	0.1ul	10	1.50%	1.50%	80102010
				5	2.50%	2.50%	
				1	4.00%	4.00%	
20001120	8-channel	5-50ul	0.5ul	50	1.00%	0.50%	80102060
				25	1.50%	1.00%	
				5	3.00%	2.00%	
20001130	8-channel	50-300ul	5ul	300	0.70%	0.25%	80102060
				150	1.00%	0.50%	
				50	1.50%	0.80%	
20001210	12-channel	0.5-10ul	0.1ul	10	1.50%	1.50%	80102010
				5	2.50%	2.50%	
				1	4.00%	4.00%	
20001220	12-channel	5-50ul	0.5ul	50	1.00%	0.50%	80102060
				25	1.50%	1.00%	
				5	3.00%	2.00%	
20001230	12-channel	50-300ul	5ul	300	0.70%	0.25%	80102060
				150	1.00%	0.50%	
				50	1.50%	0.80%	

Order Details of Multi-channel Pipette

Product Code	Channel	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
20001110	8-channel	Multi-channel Pipette	0.5-10ul	1	pcs	335	50pcs/carton	15	0.56*0.36*0.28
20001120	8-channel	Multi-channel Pipette	5-50ul	1	pcs	335	50pcs/carton	16	0.56*0.36*0.28
20001130	8-channel	Multi-channel Pipette	50-300ul	1	pcs	335	50pcs/carton	16.5	0.56*0.36*0.28
20001210	12-channel	Multi-channel Pipette	0.5-10ul	1	pcs	400	50pcs/carton	16	0.56*0.36*0.28
20001220	12-channel	Multi-channel Pipette	5-50ul	1	pcs	400	50pcs/carton	17	0.56*0.36*0.28
20001230	12-channel	Multi-channel Pipette	50-300ul	1	pcs	400	50pcs/carton	18	0.56*0.36*0.28

Pipetfocus Whole Autoclaved Single Channel Adjustable Pipette



Availability for autoclave of the whole device
Volume ranging from 0.2µl to 10ml
Easy-to-clean and replaceable parts

- A rotary button cap to ensure smooth and stable liquid transfer
- Each model has color mark with different color to show volume range
- The handle is made of high-quality materials, with an adequate thickness to minimize the impact of palm temperature on the volume.
- The tip pusher is detachable to facilitate cleaning and maintenance
- The tip is made of pure polypropylene material. No liquid droplets remain on the inner wall of the tip after expulsion. All technical performances remain unchanged after high temperature disinfection of the tip.



The Parameters of Pipetfocus								
Product Code	Model	Channel	Volume Range	Increment	Test Volume	Accuracy	Precision	Tips Code
40001010	P2	1	0.2-2ul	0.002ul	0.2ul 0.5ul 2ul	±12.00% ±5.00% ±1.50%	≤6.00% ≤2.50% ≤0.70%	80102010
40001020	P10	1	0.5-10ul	0.02ul	1ul 5ul 10ul	±2.50% ±2.00% ±1.00%	≤1.50% ≤0.80% ≤0.60%	80102010
40001030	P20	1	2-20ul	0.02ul	2ul 10ul 20ul	±5.00% ±1.00% ±1.00%	≤1.50% ≤0.60% ≤0.50%	80102060
40001040	P100	1	10-100ul	0.2ul	10ul 50ul 100ul	±1.20% ±0.80% ±0.80%	≤0.60% ≤0.30% ≤0.20%	80102060
40001050	P200	1	20-200ul	0.2ul	20ul 100ul 200ul	±2.00% ±1.00% ±1.00%	≤0.50% ≤0.40% ≤0.30%	80102060
40001060	P1000	1	100-1000ul	2ul	100ul 500ul 1000ul	±1.50% ±0.80% ±0.80%	≤0.40% ≤0.30% ≤0.20%	80102090
40001070	P5000	1	1000-5000ul	20ul	1000ul 2500ul 5000ul	±1.00% ±0.60% ±0.60%	≤0.40% ≤0.30% ≤0.20%	80102100
40001080	P10ml	1	1-10ml	20ul	1ml 5ml 10ml	±3.00% ±0.60% ±0.60%	≤0.60% ≤0.20% ≤0.16%	80102110

Order Details of Pipetfocus								
Product Code	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
40001010	Pipetfocus	0.2-2ul	1	pcs	97.00	50pcs/carton	13	0.49*0.39*0.29
40001020	Pipetfocus	0.5-10ul	1	pcs	97.00	50pcs/carton	13	0.49*0.39*0.29
40001030	Pipetfocus	2-20ul	1	pcs	97.00	50pcs/carton	13	0.49*0.39*0.29
40001040	Pipetfocus	10-100ul	1	pcs	97.00	50pcs/carton	13	0.49*0.39*0.29
40001050	Pipetfocus	20-200ul	1	pcs	97.00	50pcs/carton	13	0.49*0.39*0.29
40001060	Pipetfocus	100-1000ul	1	pcs	97.00	50pcs/carton	14	0.49*0.39*0.29
40001070	Pipetfocus	1000-5000ul	1	pcs	97.00	50pcs/carton	14.4	0.49*0.39*0.29
40001080	Pipetfocus	1-10ml	1	pcs	97.00	50pcs/carton	19.9	0.49*0.39*0.29
40002010	Pipetfocus	0.5ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002020	Pipetfocus	1ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002030	Pipetfocus	5ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002040	Pipetfocus	10ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002050	Pipetfocus	20ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002060	Pipetfocus	25ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002070	Pipetfocus	50ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002080	Pipetfocus	100ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002090	Pipetfocus	200ul	1	pcs	75.00	50pcs/carton	13	0.49*0.39*0.29
40002100	Pipetfocus	250ul	1	pcs	75.00	50pcs/carton	14	0.49*0.39*0.29
40002110	Pipetfocus	500ul	1	pcs	75.00	50pcs/carton	14	0.49*0.39*0.29
40002120	Pipetfocus	1000ul	1	pcs	75.00	50pcs/carton	14	0.49*0.39*0.29
40002130	Pipetfocus	2000ul	1	pcs	75.00	50pcs/carton	14.4	0.49*0.39*0.29
40002140	Pipetfocus	2500ul	1	pcs	75.00	50pcs/carton	14.4	0.49*0.39*0.29
40002150	Pipetfocus	5000ul	1	pcs	75.00	50pcs/carton	14.4	0.49*0.39*0.29
40002160	Pipetfocus	10ml	1	pcs	75.00	50pcs/carton	19.9	0.49*0.39*0.29

Pipette with Switch



Pipette with switch is available to assembly with an on/off switch and cable.

This can be used with coagulation analyzers or any other instruments that require precise timing.



The Parameters of Pipette with Switch

Product Code	Channel	Volume Range	Increment	Test Volume	Accuracy	Precision	Tips Code
20003300	1	0.5-10ul	0.1ul	10ul	±1.00%	±0.80%	80102010
				5ul	±1.50%	±1.50%	
				1ul	±2.50%	±1.50%	
20003310	1	2-20ul	0.5ul	20ul	±0.90%	±0.40%	80102060
				10ul	±1.20%	±1.00%	
				2ul	±3.00%	±2.00%	
20003320	1	5-50ul	0.5ul	50ul	±0.60%	±0.30%	80102060
				25ul	±0.90%	±0.60%	
				5ul	±2.00%	±2.00%	
20003330	1	10-100ul	1ul	100ul	±0.80%	±0.15%	80102060
				50ul	±1.00%	±0.40%	
				10ul	±3.00%	±1.50%	
20003340	1	20-200ul	1ul	200ul	±0.60%	±0.15%	80102060
				100ul	±0.80%	±0.30%	
				20ul	±3.00%	±1.00%	
20003350	1	100-1000ul	5ul	1000ul	±0.60%	±0.20%	80102090
				500ul	±0.70%	±0.25%	
				100ul	±1.50%	±0.70%	
20003360	1	1000-5000ul	50ul	5000ul	±0.50%	±0.15%	80102100
				2500ul	±0.60%	±0.30%	
				1000ul	±0.70%	±0.30%	

Order Details of Pipette with Switch

Product Code	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
20003300	Pipette with Switch	0.5-10ul	1	pcs	88.00	50pcs/carton	10.5	0.475*0.39*0.28
20003310	Pipette with Switch	2-20ul	1	pcs	88.00	50pcs/carton	10.5	0.475*0.39*0.28
20003320	Pipette with Switch	5-50ul	1	pcs	88.00	50pcs/carton	10.5	0.475*0.39*0.28
20003330	Pipette with Switch	10-100ul	1	pcs	88.00	50pcs/carton	10.7	0.475*0.39*0.28
20003340	Pipette with Switch	20-200ul	1	pcs	88.00	50pcs/carton	10.75	0.475*0.39*0.28
20003350	Pipette with Switch	100-1000ul	1	pcs	88.00	50pcs/carton	11.2	0.475*0.39*0.28
20003360	Pipette with Switch	1-5ml	1	pcs	88.00	50pcs/carton	11.3	0.475*0.39*0.28
20003410	Fixed Pipette with Switch	0.5ul	1	pcs	70.00	50pcs/carton	10.5	0.475*0.39*0.28
20003420	Fixed Pipette with Switch	1ul	1	pcs	70.00	50pcs/carton	10.5	0.475*0.39*0.28
20003430	Fixed Pipette with Switch	5ul	1	pcs	70.00	50pcs/carton	10.5	0.475*0.39*0.28
20003440	Fixed Pipette with Switch	10ul	1	pcs	70.00	50pcs/carton	10.5	0.475*0.39*0.28
20003450	Fixed Pipette with Switch	20ul	1	pcs	70.00	50pcs/carton	10.5	0.475*0.39*0.28
20003460	Fixed Pipette with Switch	25ul	1	pcs	70.00	50pcs/carton	10.5	0.475*0.39*0.28
20003470	Fixed Pipette with Switch	50ul	1	pcs	70.00	50pcs/carton	10.5	0.475*0.39*0.28
20003480	Fixed Pipette with Switch	100ul	1	pcs	70.00	50pcs/carton	10.75	0.475*0.39*0.28
20003490	Fixed Pipette with Switch	200ul	1	pcs	70.00	50pcs/carton	10.75	0.475*0.39*0.28
20003500	Fixed Pipette with Switch	250ul	1	pcs	70.00	50pcs/carton	11.2	0.475*0.39*0.28
20003510	Fixed Pipette with Switch	500ul	1	pcs	70.00	50pcs/carton	11.2	0.475*0.39*0.28
20003520	Fixed Pipette with Switch	1000ul	1	pcs	70.00	50pcs/carton	11.2	0.475*0.39*0.28
20003530	Fixed Pipette with Switch	2000ul	1	pcs	70.00	50pcs/carton	11.3	0.475*0.39*0.28
20003540	Fixed Pipette with Switch	2500ul	1	pcs	70.00	50pcs/carton	11.3	0.475*0.39*0.28
20003550	Fixed Pipette with Switch	5000ul	1	pcs	70.00	50pcs/carton	11.3	0.475*0.39*0.28

Coloured Adjustable Pipette-Five Fixed Volume



The 5 colors above indicate different volume ranges respectively
Volume range can be custom-made

- A left suction button with volume indicator
- The volumes for adjustment are shown on the panel .
- Availability for self calibration of volume through the window
- An identification sheath with different color, for clearer recognition of different specifications
- A tip pusher made of stainless steel material, with excellent chemical corrosion resistance
- The tip is made of pure polypropylene material to withstand a high temperature of 121°C. No liquid droplets remain on the inner wall of the tip after expulsion.

Volume range can be custom-made



The Parameters of Colored Adjustable & Fixed Pipettes

Volume	Accuracy	Precision	Volume	Accuracy	Precision	Volume	Accuracy	Precision
5ul	±8.00%	≤5.00%	40ul	±3.00%	≤1.50%	300ul	±2.00%	±1.00%
10ul	±4.80%	≤2.50%	50ul	±3.00%	≤1.50%	400ul	±2.00%	±1.00%
15ul	±4.80%	≤2.50%	100ul	±2.00%	≤1.00%	500ul	±2.00%	±1.00%
20ul	±3.00%	≤1.50%	150ul	±2.00%	≤1.00%	600ul	±2.00%	±1.00%
25ul	±3.00%	≤1.50%	200ul	±2.00%	≤1.00%	800ul	±1.50%	≤0.80%
30ul	±3.00%	≤1.50%	250ul	±2.00%	≤1.00%	1000ul	±1.50%	≤0.80%

The Specifications and Order Details of Coloured Adjustable Pipette

Product Code	Channel	Volume Range	Increment	Qty	Unit	Retail Price (EURO)	Tips Code	Package	Gross weight (kg)	Dimensions in cubic meters
10002330	1	5-25ul	5ul	1	pcs	35	80102060	50pcs/carton	8.95	0.475*0.39*0.28
10002340	1	10-50ul	10ul	1	pcs	35	80102060	50pcs/carton	9.1	0.475*0.39*0.28
10002350	1	50-250ul	50ul	1	pcs	35	80102060	50pcs/carton	9.6	0.475*0.39*0.28
10002360	1	100-500ul	100ul	1	pcs	35	80102090	50pcs/carton	9.7	0.475*0.39*0.28
10002370	1	200-1000ul	200ul	1	pcs	35	80102090	50pcs/carton	10.1	0.475*0.39*0.28
10002230	1	5ul		1	pcs	30	80102060	50pcs/carton	8.95	0.475*0.39*0.28
10002240	1	10ul		1	pcs	30	80102060	50pcs/carton	8.95	0.475*0.39*0.28
10002250	1	20ul		1	pcs	30	80102060	50pcs/carton	8.95	0.475*0.39*0.28
10002260	1	25ul		1	pcs	30	80102060	50pcs/carton	8.95	0.475*0.39*0.28
10002270	1	50ul		1	pcs	30	80102060	50pcs/carton	9.4	0.475*0.39*0.28
10002280	1	100ul		1	pcs	30	80102060	50pcs/carton	9.4	0.475*0.39*0.28
10002290	1	200ul		1	pcs	30	80102060	50pcs/carton	9.4	0.475*0.39*0.28
10002300	1	250ul		1	pcs	30	80102060	50pcs/carton	9.6	0.475*0.39*0.28
10002310	1	500ul		1	pcs	30	80102090	50pcs/carton	9.7	0.475*0.39*0.28
10002320	1	1000ul		1	pcs	30	80102090	50pcs/carton	10.1	0.475*0.39*0.28

Purple Adjustable Pipette-Five Fixed Volume



Five fixed volume ranges with convenience of operation
Switching from the minimum volume range to the maximum volume range can be carried out within two seconds

- Five clear, fixed volume scales
- A deft handle design
- A tip pusher made of stainless steel material, with excellent chemical corrosion resistance



The Parameters of Purple Adjustable & Fixed Pipette

Volume	Accuracy	Precision	Volume	Accuracy	Precision	Volume	Accuracy	Precision
5ul	±8.00%	≤5.00%	100ul	±2.00%	≤1.00%	800ul	±1.50%	≤0.80%
10ul	±4.80%	≤2.50%	150ul	±2.00%	≤1.00%	1000ul	±1.50%	≤0.80%
15ul	±4.80%	≤2.50%	200ul	±2.00%	≤1.00%	1500ul	±1.00%	≤0.50%
20ul	±3.00%	≤1.50%	250ul	±2.00%	≤1.00%	2000ul	±1.00%	≤0.50%
25ul	±3.00%	≤1.50%	300ul	±2.00%	±1.00%	2500ul	±1.00%	≤0.50%
30ul	±3.00%	≤1.50%	400ul	±2.00%	±1.00%	3000ul	±1.00%	≤0.50%
40ul	±3.00%	≤1.50%	500ul	±2.00%	±1.00%	4000ul	±1.00%	≤0.50%
50ul	±3.00%	≤1.50%	600ul	±2.00%	±1.00%	5000ul	±1.00%	≤0.50%

Order Details of Purple Adjustable & Fixed Pipette

Product Code	Channel	Volume Range	Increment	Qty	Unit	Retail Price (EURO)	Tips Code	Package	Gross weight (kg)	Dimensions in cubic meters
10001130	1	5-25ul	5ul	1	pcs	35	80102030	50pcs/carton	6.05	0.475*0.39*0.28
10001140	1	10-50ul	10ul	1	pcs	35	80102030	50pcs/carton	6.1	0.475*0.39*0.28
10001150	1	50-250ul	50ul	1	pcs	35	80102030	50pcs/carton	6.3	0.475*0.39*0.28
10001160	1	100-500ul	100ul	1	pcs	35	80102080	50pcs/carton	6.7	0.475*0.39*0.28
10001170	1	200-1000ul	200ul	1	pcs	35	80102080	50pcs/carton	6.95	0.475*0.39*0.28
10001180	1	500-2500ul	500ul	1	pcs	37	80102100	50pcs/carton	7.4	0.475*0.39*0.28
10001185	1	1000-5000ul	1000ul	1	pcs	37	80102100	50pcs/carton	8.1	0.475*0.39*0.28
10001010	1	5ul		1	pcs	30	80102030	50pcs/carton	6.15	0.475*0.39*0.28
10001020	1	10ul		1	pcs	30	80102030	50pcs/carton	6.15	0.475*0.39*0.28
10001030	1	20ul		1	pcs	30	80102030	50pcs/carton	6.2	0.475*0.39*0.28
10001040	1	25ul		1	pcs	30	80102030	50pcs/carton	6.2	0.475*0.39*0.28
10001050	1	50ul		1	pcs	30	80102030	50pcs/carton	6.35	0.475*0.39*0.28
10001060	1	100ul		1	pcs	30	80102030	50pcs/carton	6.35	0.475*0.39*0.28
10001070	1	200ul		1	pcs	30	80102030	50pcs/carton	6.35	0.475*0.39*0.28
10001080	1	250ul		1	pcs	30	80102030	50pcs/carton	6.35	0.475*0.39*0.28
10001090	1	500ul		1	pcs	30	80102080	50pcs/carton	7	0.475*0.39*0.28
10001100	1	1000ul		1	pcs	30	80102080	50pcs/carton	7	0.475*0.39*0.28
10001110	1	2000ul		1	pcs	32	80102100	50pcs/carton	7.95	0.475*0.39*0.28
10001120	1	2500ul		1	pcs	32	80102100	50pcs/carton	8.1	0.475*0.39*0.28
10001125	1	5000ul		1	pcs	32	80102100	50pcs/carton	8.3	0.475*0.39*0.28

Dispensing Pipette



It features a light weight design with ease of use. One suction allows for continuously dispensing 45 times. It has strong universality and is able to adapt to 7 types of dispensing pipette tip.

The dispensing range is 10-5000µl and for each type of tip, five dispensing volumes are available for selection. Only attach a suitable tip and set the dispensing volume required, and then enjoy the fun of fast dispensing. A unique ergonomic handle design brings comfort of use. The light weight makes it suitable for one-hand operation.

Rongtai Dispensing Pipette employs the external piston replacement principle and is especially suitable for dispensing corrosive and sticky liquids.

Tips for dispensing pipette are packed separately. For use of a 25 or 50ml dispensing pipette tip, an adaptor is required. Each dispensing pipette is provided with an adaptor in the package.



The Parameters of Dispensing Pipette

Tips Code	Tips Volume	Test Volume	Accuracy		Precision	
90101010	0.5ml	50ul	±0.5ul	±1.0%	0.3ul	0.6%
		10ul	±0.1ul	±1.0%	0.1ul	1.0%
90101020	1.25ml	125ul	±1.0ul	±0.8%	0.6ul	0.5%
		25ul	±0.2ul	±0.8%	0.2ul	0.8%
90101030	2.5ml	250ul	±2.0ul	±0.8%	0.8ul	0.3%
		50ul	±0.4ul	±0.8%	0.4ul	0.7%
90101040	5.0ml	500ul	±4.0ul	±0.8%	1.0ul	0.2%
		100ul	±0.8ul	±0.8%	0.5ul	0.5%
90101050	12.5ml	1250ul	±7.5ul	±0.6%	2.5ul	0.2%
		250ul	±1.5ul	±0.6%	1.0ul	0.4%
90101060	25ml	2500ul	±12.5ul	±0.5%	7.5ul	0.3%
		500ul	±2.5ul	±0.5%	2.5ul	0.5%
90101070	50ml	5000ul	±25.0ul	±0.5%	15.0ul	0.3%
		1000ul	±5.0ul	±0.5%	5.0ul	0.5%

Order Details of Dispensing Pipette

Product Code	Volume Range	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
50001005	10-5000ul	1	pcs	200.00	50pcs/carton	10	0.475*0.39*0.28

Order Details of Tips for Dispensing Pipette

Product Code	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
90101010	Tips for Dispensing Pipette	0.5ml	1	pcs	3.65	100pcs/box	0.25	0.22*0.135*0.115
90101020	Tips for Dispensing Pipette	1.25ml	1	pcs	3.65	100pcs/box	0.25	0.22*0.135*0.115
90101030	Tips for Dispensing Pipette	2.5ml	1	pcs	3.65	100pcs/box	0.25	0.22*0.135*0.115
90101040	Tips for Dispensing Pipette	5.0ml	1	pcs	3.15	50pcs/box	0.25	0.22*0.135*0.115
90101050	Tips for Dispensing Pipette	12.5ml	1	pcs	3.15	50pcs/box	0.25	0.22*0.135*0.115
90101060	Tips for Dispensing Pipette	25ml	1	pcs	3.95	20pcs/box	0.25	0.22*0.135*0.115
90101070	Tips for Dispensing Pipette	50ml	1	pcs	4.25	10pcs/box	0.25	0.22*0.135*0.115

Adjustable Glass-Injection Dispenser

Availability for autoclave of the whole device
A non-fragile, durable silicone tube made of special materials
Convenience of use with high efficiency, and low repetitive error, capable of greatly reducing time and costs in research and experiment

- Solvent Resistance: general acid and alkali
- Material: PP



The Parameters of Adjustable Glass-Injection Dispenser

Product Code	Description	Volume	Test Volume	Accuracy
30011420	Adjustable Glass-Injection Dispenser	0.38ml	0.38ml	±0.01
30011430	Adjustable Glass-Injection Dispenser	0.1-1ml	0.1ml	±0.02
30011440	Adjustable Glass-Injection Dispenser	0.5-5ml	0.5ml	±0.04
30011450	Adjustable Glass-Injection Dispenser amber glass	0.5-5ml	0.5ml	±0.04
30011460	Adjustable Glass-Injection Dispenser translucent glass	1-10ml	1ml	±0.1
30011470	Adjustable Glass-Injection Dispenser amber glass	1-10ml	1ml	±0.1
30011480	Adjustable Glass-Injection Dispenser plastic	1-10ml	1ml	±0.1

Order Details of Single Channel Digital Variable Pipette

Product Code	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
30011420	Adjustable Glass-Injection Dispenser	0.38ml	1	pcs	14	10pcs/carton	3.7	0.52*0.28*0.29
30011430	Adjustable Glass-Injection Dispenser	0.1-1ml	1	pcs	14	10pcs/carton	3.9	0.52*0.28*0.29
30011440	Adjustable Glass-Injection Dispenser	0.5-5ml	1	pcs	14	10pcs/carton	3.9	0.52*0.28*0.29
30011450	Adjustable Glass-Injection Dispenser amber glass	0.5-5ml	1	pcs	15	10pcs/carton	4.2	0.62*0.25*0.315
30011460	Adjustable Glass-Injection Dispenser translucent glass	1-10ml	1	pcs	18	10pcs/carton	6.5	0.66*0.315*0.335
30011470	Adjustable Glass-Injection Dispenser amber glass	1-10ml	1	pcs	18	10pcs/carton	6	0.66*0.315*0.335
30011480	Adjustable Glass-Injection Dispenser plastic	1-10ml	1	pcs	18	10pcs/carton	2.9	0.66*0.315*0.335

Bottle-Top Dispenser

5-25ml Built-In Bottle-Top Dispenser is suitable for continuously dispensing liquid with a fixed volume in laboratories, and is applicable in petrochemical industry, medical and health care, metallurgy, scientific research, environmental protection, and other fields.

- Solvent Resistance: general acid and alkali
- Material: ABS



The Parameters of Bottle-Top Dispenser

Product Code	Description	Volume	Test Volume	Accuracy
30011490	Bottle-Top Dispenser	5-25ml	5ml	±0.7

Order Details of Bottle-Top Dispenser

Product Code	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
30011490	Bottle-Top Dispenser	5-25ml	1	pcs	32	10pcs/carton	6.1	0.72*0.21*0.2

Pen PH Recording



The PHB-1 meter uses integrated circuit to design, it is used to measure pH of liquids. The meter widely applies to environmental monitoring, swimming pool, food processing and other industrial applications.

SPECIFICATIONS:

Model: PHB-1
 pH Range: 0.0~14.0 pH
 Resolution: 0.1 pH
 Operation Temperature: 5 to 35°C
 Power: 3V batteries × 3
 Dimensions: 130 × 28 × 18 mm
 Weight: 100g
 Package: 50pcs/carton
 Product Code: 80208330

Retail Price(EURO): 60.00/pcs

Microliter Syringe



Order Details of Microliter Syringe

Product Code	Description	Volume	Qty	Unit	Retail Price (EURO)	Package	Gross weight (kg)	Dimensions in cubic meters
30012005	MICROLITER SYRINGE GC	0.5ul	1	pcs	15	100pcs/carton	6	0.46*0.37*0.19
30012010	MICROLITER SYRINGE GC	1ul	1	pcs	15	100pcs/carton	6	0.46*0.37*0.19
30012015	MICROLITER SYRINGE GC	5ul	1	pcs	15	100pcs/carton	6	0.46*0.37*0.19
30012020	MICROLITER SYRINGE GC	10ul	1	pcs	13	100pcs/carton	6	0.46*0.37*0.19
30012030	MICROLITER SYRINGE GC	25ul	1	pcs	13	100pcs/carton	6.1	0.46*0.37*0.19
30012035	MICROLITER SYRINGE GC	50ul	1	pcs	13	100pcs/carton	6.1	0.46*0.37*0.19
30012040	MICROLITER SYRINGE GC	100ul	1	pcs	13	100pcs/carton	6.2	0.46*0.37*0.19
30012045	MICROLITER SYRINGE GC	250ul	1	pcs	30	100pcs/carton	7	0.46*0.37*0.19
30012050	MICROLITER SYRINGE GC	500ul	1	pcs	30	100pcs/carton	7.3	0.46*0.37*0.19
30012055	MICROLITER SYRINGE GC	1000ul	1	pcs	30	100pcs/carton	7.3	0.46*0.37*0.19
30012006	MICROLITER SYRINGE LC	0.5ul	1	pcs	20	100pcs/carton	6	0.46*0.37*0.19
30012011	MICROLITER SYRINGE LC	1ul	1	pcs	20	100pcs/carton	6	0.46*0.37*0.19
30012016	MICROLITER SYRINGE LC	5ul	1	pcs	20	100pcs/carton	6	0.46*0.37*0.19
30012021	MICROLITER SYRINGE LC	10ul	1	pcs	19	100pcs/carton	6	0.46*0.37*0.19
30012031	MICROLITER SYRINGE LC	25ul	1	pcs	19	100pcs/carton	6.1	0.46*0.37*0.19
30012036	MICROLITER SYRINGE LC	50ul	1	pcs	19	100pcs/carton	6.1	0.46*0.37*0.19
30012041	MICROLITER SYRINGE LC	100ul	1	pcs	19	100pcs/carton	6.2	0.46*0.37*0.19
30012046	MICROLITER SYRINGE LC	250ul	1	pcs	33	100pcs/carton	7	0.46*0.37*0.19
30012051	MICROLITER SYRINGE LC	500ul	1	pcs	33	100pcs/carton	7.3	0.46*0.37*0.19
30012056	MICROLITER SYRINGE LC	1000ul	1	pcs	33	100pcs/carton	7.3	0.46*0.37*0.19

Enzyme Label Plate



- Opaque walls to prevent well-to-well crosstalk
- Optically clear flat well bottom permits direct microscopic viewing
- Can be used for both top and bottom reading instruments
- Flat bottoms with 250µL total volume
- Recommended working volumes of 50 to 150µL
- Treated for optimal cell attachment
- Individual alphanumeric codes for well identification
- Bulk packed 100 per box

Order Details of Enzyme Label Plate

Product Code	Description	Specification	Qty	Unit	Retail Price	Package	Gross Weight (KG)	Dimensions
80208109	Enzyme Label Plate (Detachable) 96-well 8 well*12	Flat Bottom/Medium Binding	1	pcs	2.55	100pcs/carton	5.4	0.375*0.27*0.18
80208115	Enzyme Label Plate (Detachable) 96-well 8 well*12	Flat Bottom/High Binding	1	pcs	3.05	100pcs/carton	5.4	0.375*0.27*0.18
80208116	Enzyme Label Plate (Detachable) 96-well 8 well*12	Flat Bottom/Aminated	1	pcs	2.95	100pcs/carton	5.4	0.375*0.27*0.18
80208110	Enzyme Label Plate (Detachable) 96-well 12 well*8	Flat Bottom/Medium Binding	1	pcs	2.85	100pcs/carton	5.4	0.375*0.27*0.18
80208117	Enzyme Label Plate (Detachable) 96-well 12 well*8	Flat Bottom/High Binding	1	pcs	3.05	100pcs/carton	5.4	0.375*0.27*0.18
80208132	Enzyme Label Plate (Detachable) 96-well 12 well*8	Flat Bottom/Aminated	1	pcs	3.99	100pcs/carton	5.4	0.375*0.27*0.18
80208131	Enzyme Label Plate (Detachable) 96-well 8 well*12 Sterilized	Flat Bottom/Medium Binding	1	pcs	3.55	100pcs/carton	5.4	0.375*0.27*0.18
80208118	Enzyme Label Plate (Detachable) 96-well 8 well*12 Sterilized	Flat Bottom/High Binding	1	pcs	2.86	100pcs/carton	5.4	0.375*0.27*0.18
80208119	Enzyme Label Plate (Detachable) 96-well 8 well*12 Sterilized	Flat Bottom/Aminated	1	pcs	2.86	100pcs/carton	5.4	0.375*0.27*0.18
80208120	Enzyme Label Plate (Detachable) 96-well 12 well*8 Sterilized	Flat Bottom/Medium Binding	1	pcs	2.86	100pcs/carton	5.4	0.375*0.27*0.18
80208121	Enzyme Label Plate (Detachable) 96-well 12 well*8 Sterilized	Flat Bottom/High Binding	1	pcs	2.96	100pcs/carton	5.4	0.375*0.27*0.18
80208122	Enzyme Label Plate (Detachable) 96-well 12 well*8 Sterilized	Flat Bottom/Aminated	1	pcs	2.86	100pcs/carton	5.4	0.375*0.27*0.18
80208111	Enzyme Label Strip	8 well/Medium Binding	1	strip	0.3	4000pcs/carton	11	0.58*0.28*0.2
80208123	Enzyme Label Strip	8 well/High Binding	1	strip	0.3	4000pcs/carton	11	0.58*0.28*0.2
80208125	Enzyme Label Strip	8 well/Aminated	1	strip	0.3	4000pcs/carton	11	0.58*0.28*0.2
80208112	Enzyme Label Strip	12 well/Medium Binding	1	strip	0.3	3000pcs/carton	11	0.58*0.28*0.2
80208124	Enzyme Label Strip	12 well/High Binding	1	strip	0.3	3000pcs/carton	11	0.58*0.28*0.2
80208126	Enzyme Label Strip	12 well/Aminated	1	strip	0.3	3000pcs/carton	11	0.58*0.28*0.2
80208113	Enzyme Label Strip Sterilized	8 well/Medium Binding	1	strip	0.35	4000pcs/carton	11	0.58*0.28*0.2
80208127	Enzyme Label Strip Sterilized	8 well/High Binding	1	strip	0.35	4000pcs/carton	11	0.58*0.28*0.2
80208129	Enzyme Label Strip Sterilized	8 well/Aminated	1	strip	0.35	4000pcs/carton	11	0.58*0.28*0.2
80208114	Enzyme Label Strip Sterilized	12 well/Medium Binding	1	strip	0.35	3000pcs/carton	11	0.58*0.28*0.2
80208128	Enzyme Label Strip Sterilized	12 well/High Binding	1	strip	0.35	3000pcs/carton	11	0.58*0.28*0.2
80208130	Enzyme Label Strip Sterilized	12 well/Aminated	1	strip	0.35	3000pcs/carton	11	0.58*0.28*0.2

Tips for Pipette



Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102010	3.5*32mm	10ul	1000pcs/bag	100bag/carton	19.5/bag



Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102030	4*49mm	200ul	500pcs/bag	100bag/carton	8.9/bag



Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102050	5.5*50mm	200ul	500pcs/bag	100bag/carton	8.9/bag



Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102060	5.5*52mm	300ul	500pcs/bag	100bag/carton	8.9/bag



Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102080	7.4*70mm	1000ul	500pcs/bag	20bag/carton	19.9/bag



Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102090	7.5*70mm	1000ul	500pcs/bag	20bag/carton	19.9/bag



Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102100	13*120mm	5ml	50pcs/bag	100bag/carton	8.9/bag




Product Code	Specification	Volume	Qty	Package	Price(EURO)
80102110	16*150mm	10ml	100pcs/bag	50bag/carton	32/bag

Racks for Tips & Centrifuge Tubes



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208010	Racks for Tips 1000ul	60 hole	1	pcs	5.5



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208011	Racks for Tips 1000ul	100 hole	1	pcs	7



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208020	Racks for Tips 300ul	96 hole	1	pcs	5.5



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208023	Racks for Tips 300ul	96 hole	1	pcs	5.5




Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208025	Racks for Tips 10ul	96 hole	1	pcs	5.5



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208030	Racks for Centrifuge Tube 1.5/2.0ml	72 hole	1	pcs	8
80208035	Racks for Centrifuge Tube 0.5ml	72 hole	1	pcs	8



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208062	Multifunctional Board for Centrifuge Tube	1.5/2.0ml	1	pcs	7



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208043	Multifunctional Rack for Centrifuge Tube	50ml / 15ml	1	pcs	9

Pipette Stander & Tube Rack



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208044	Pipettor Solution Basins	50ml	1	pcs	1.5



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208050	Pipette Stander	Detachable	1	pcs	20




Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208045	Pipette Stander	Acryl	1	pcs	29



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208055	Pipette Stander	Multifunction	1	pcs	29



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208056	Pipette Stander	for Digital Pipette	1	pcs	40



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208037	Stainless Steel Rack for Test Tube	ø13*50 hole	1	pcs	18
80208038	Stainless Steel Rack for Test Tube	ø16*50 hole	1	pcs	18
80208039	Stainless Steel Rack for Test Tube	ø18*50 hole	1	pcs	18



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208060	Multifunctional Rack for Test Tube	50 hole	1	pcs	12



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208040	Plastic Rack for Test Tube	ø13*50 hole	1	pcs	7
80208041	Plastic Rack for Test Tube	ø16*50 hole	1	pcs	7
80208042	Plastic Rack for Test Tube	ø18*50 hole	1	pcs	7

Centrifuge Tubes



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100020	Centrifuge Tubes	1.5ml	500pcs	bag	12.5
50100025	Centrifuge Tubes	2ml	500pcs	bag	12.5



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100140	Centrifuge Tubes with Screw Cap	100ml	50pcs	bag	33
50100150	Centrifuge Tubes with Cap	50ml	50pcs	bag	15
50100060	Centrifuge Tubes with Cap	100ml	50pcs	bag	25



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100040	Centrifuge Tubes	10ml	200pcs	bag	20



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100130	Centrifuge Tubes with Screw Cap	50ml	50pcs	bag	12



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100220	Cryovial Tubes	5ml	200pcs	bag	33



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100210	Cryovial Tubes	1.8ml	500pcs	bag	42



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100005	Centrifuge Tubes	0.2ml	1000pcs	bag	17
50100010	Centrifuge Tubes	0.5ml	1000pcs	bag	17



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100030	Centrifuge Tubes	5ml	300pcs	bag	15



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
50100120	Centrifuge Tubes with Screw Cap	15ml	100pcs	bag	20

Others



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208224	Cell Culture Plate	6 hole	1	pcs	2.9
80208239	Cell Culture Plate	12 hole	1	pcs	2.9
80208225	Cell Culture Plate	24 hole	1	pcs	2.9
80208226	Cell Culture Plate	96 hole	1	pcs	2.9



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208227	Petri Dishes	ø35	1	pcs	0.5
80208228	Petri Dishes	ø60	1	pcs	0.5
80208229	Petri Dishes	ø70	1	pcs	0.5
80208230	Petri Dishes	ø90	1	pcs	0.5
80208237	Petri Dishes	ø100	1	pcs	0.85
80208238	Petri Dishes	ø120	1	pcs	0.95
80208231	Petri Dishes	ø150	1	pcs	0.95



Product Code	Description	Qty	Unit	Retail Price (EURO)
80208160	Counting Chambers	1	pcs	8



Product Code	Description	Qty	Unit	Retail Price (EURO)
80208150	Implant Inoculable Needle	1	pcs	3



Product Code	Description	Qty	Unit	Retail Price (EURO)
80208320	Nickel Chrome Wire	20 wires	pack	1.2



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208091	Hemagglutination Board(Acryl)	V110°	1	pcs	4.5
80208092	Hemagglutination Board(Acryl)	V90°	1	pcs	4.5
80208093	Hemagglutination Board(Acryl)	U	1	pcs	4.5



Product Code	Description	Specification	Qty	Unit	Retail Price (EURO)
80208101	Hemagglutination Board(PVC)	V110°	1	pcs	0.4
80208100	Hemagglutination Board(PVC)	U	1	pcs	0.4